

Experimental methods I

Example experimental set ups and methods

M.P. Vaughan




Example experimental setup (1)



- Set up to investigate the dynamics of laser pairs (mounted on the same chip) subject to mutual optical injection
- Acronyms used:
 - DBR – Distributed Bragg Reflector (laser)
 - TEC – Thermo-Electric Controller
 - FPI – Fabry Perot Interferometer
 - ESA – Electrical Spectrum Analyser
 - OSA – Optical Spectrum Analyser

Example experimental setup (1)




Schematic of the experimental setup for the mutual injection measurements.

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Laser mounting



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Fibre coupling to semiconductor laser



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Elements of experimental set up




- Positioners
 - Translation stages (x, y & z)
 - Rotation stages
 - Tilt stages

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Nano-stage positioner




Micrometer screws
with Vernier scales

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Vernier scale




Vernier scale

Photograph of a Vernier caliper. Copyright ArtMechanic
(2004)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Vernier scale – close up




Close-up of a Vernier scale. Copyright Glenn (2005)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Micrometer – metric scale



Photograph of a micrometer. Copyright Glenn McKechnie (2005)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Micrometer – Vernier scale



Photograph of a micrometer. Copyright Glenn McKechnie (2005)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Elements of experimental set up




- Laser peripherals
 - Current source
 - Temperature controller (in conjunction with a thermistor)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Laser peripherals



Thermistor
(resistance varies
with temperature)

Thermo-
Electric
Controller

N.B. Thermistor generally needs to be calibrated

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Elements of experimental set up




- Fibre optics
 - Optical fibre
 - Lens-ended fibre (for coupling laser light)
 - Optical isolator (an optical 'diode')
 - Polarisation controllers ('bunny ears')
 - Couplers (coupling / splitting light from fibre)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Optical isolator



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Optical isolator - Faraday rotation



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Optical isolator - Faraday rotation



T. Amemiya and Y. Nakano (2010). Single Mode Operation of 1.5- μm Waveguide Optical Isolators Based on the Nonreciprocal-loss Phenomenon, *Advances in Optical and Photonic Devices*, Ki Young Kim (Ed.), ISBN: 978-953-7619-76-3, InTech, DOI: 10.5772/7143. Available from: <http://www.intechopen.com/books/advances-in-optical-and-photonic-devices/single-mode-operation-of-1-5-micro-m-waveguide-optical-isolators-based-on-the-nonreciprocal-loss-phe>

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Polarisation controller ('bunny ears')




Twisting the fibre creates strain-induced **birefringence** in the fibre (light polarised in different directions travels at different speeds). The lengths of loops on the 'bunny ears' are chosen to reproduce the effect of **wave plates**.

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Fibre couplers



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Elements of experimental set up




- Optical analysis
 - Scanning Fabry Perot Interferometer
 - Optical Spectrum Analyser
 - Electrical Spectrum Analyser
 - In conjunction with a fast photodiode
 - Often termed an RF (radio frequency) analyser

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Scanning Fabry-Perot Interferometer




Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Scanning Fabry-Perot Interferometer



The FPI is a resonant cavity that may be used to measure the line-width of a laser.



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Scanning Fabry-Perot Interferometer




$\Delta\lambda$ is the **full spectral range**

$\delta\lambda$ is the **full width half maximum** of the FPI

F is the **finesse** defined by

$$F = \frac{\Delta\lambda}{\delta\lambda}.$$



Note that $\delta\lambda$ must be less than the line-width of the laser to be measured.

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Optical spectrum analyser



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Optical spectrum analyser




- Usually employs a **diffraction grating** to resolve different wavelengths
- Not as accurate as FPI for line-width
- Better for measuring wavelength / frequency over design ranges (usually optical to IR)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Electrical (RF) spectrum analyser



Performs a Fourier analysis of the electrical input signal

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Example experimental methods



- Measuring laser line-width
 - Heterodyning
 - Self-heterodyning


Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Heterodyning set up



2 laser heterodyning set up



N.B. ECL – External Cavity Laser


Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Self heterodyning set up




Self-heterodyning set up



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Line-width of ECL via self-heterodyning



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Line-width of ECL via self-heterodyning



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Example experimental setup (2)



- Investigation of novel THz devices based on photomixing
- Uses some of the same equipment for the mutual optical injection measurements
- Further acronyms:
 - DFB – Distributed FeedBack (laser)
 - EDFA – Erbium Doped Fibre Amplifier (optical amplifier)

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Photomixing




Mixing of two de-tuned optical sources

$$E = E_1 e^{i\omega_1 t} + E_2 e^{i\omega_2 t + \phi}$$

where

$$\omega_1 - \omega_2 = \delta\omega$$



In a square law detector, response is proportional to the intensity:


$$|E|^2 = (E_1^2 + E_2^2) \left(1 + \frac{2E_1 E_2}{E_1^2 + E_2^2} \cos(\delta\omega t - \phi) \right)$$

Seek beat frequency in THz

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Example experimental setup (2)



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Lens-ended fibre and THz device



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Golay cell




Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Chopper




Used to modulate optical input for lock-in amplifier



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Lock-in amplifier (and scope)



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Lock-in amplifier basics




- Allows detection of weak signals buried in noise
 - Input signal multiplied by a reference signal modulated at the same frequency
 - Leads to DC component proportional to input signal amplitude
 - Can be optimised by adjusting the phase
 - Other frequencies filtered out using band-pass filters

Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

THz power



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

Example experimental setup (3)




- Set up to measure detector quantum efficiency using correlated photons
- Elements of set up controlled centrally using LabView


Example experimental setup (3)



Basic concept




Nonlinear crystal – parametric down conversion



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Example experimental setup (3)



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Light-tight environment



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Nonlinear crystal



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Trigger photodetector (APD)



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


Device under test (DUT) in box



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics


LabView – aligning and scanning



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

LabView – coincidence counts



Coláiste na hOllscoile Corcaigh, Éire
University College Cork, Ireland

ROINN NA FISICE
Department of Physics

References to examples



M.P. Vaughan *et al*,
Mutual optical injection in coupled DBR laser pairs
Opt. Express **17**, 2033 (2009)

ID Henning, Y Sun, MP Vaughan *et al*,
Widely tuneable optoelectronic source of continuous-wave terahertz radiation
Electron. Lett. **45**, 1252 (2009) (Journal)

ID Henning, MJ Adams, MP Vaughan *et al*,
Novel antenna-integrated photodiodes with strained absorbers designed for use as terahertz sources
IEEE J. Sel. Top. Quant. **17**, 202 (2010)

J.Y. Cheung, M.P. Vaughan, J.R. Mountford and C.J. Chunnillall,
Correlated photon metrology of detectors and sources
Proc. Spie **5161**, 365-376 (2004)